



ANALYSIS

We offer a comprehensive range of testing methods for sweeteners and food supplements in tablet, powder or liquid form.

Our laboratory system documents all the obtained results and complies to country-specific requirements. Our certificates are established in German, English or French. Our experienced and knowledgeable analysts validated the methods in compliance with ICH guidelines.

SERVICES

Determination of content

- HPLC (Detection: UV, DAD)
- Automated titrator (potentiometric electrode and phototrode)

Identity

- FTIR with ATR

Physical / chemical parameters

- Measurement of loss on drying (drying chamber, halogen dryer)
 - Breaking strength and friability measurement
 - Particle size distribution (by sieve)
 - Water activity
-

